

## HOT SALE

### 980nm+1470nm Diode laser system


**10** Functions in 1 Laser







# Tırnak Mantarı Temizleme Fonksiyonu



## Lipolysis Handpiece



### Technical Data:

- Include 5 pieces cannula
- Suitable for 200~600um fiber





# Endovenous Laser Treatment (EVLT)

Colemed  
pharma

Radyal fiber ile EVLT işlevi

Using a laser for varicose veins eliminates much of the bruising and pain often associated with other vein surgeries.


HERE'S HOW IT WORKS:

-  1 A local anesthetic is applied
-  2 A small incision is made
-  3 The fiber, which emits laser energy, is inserted
-  4 The fiber is pulled back through the vein, which causes the damaged vein to collapse and seal shut
- 5 YOU'LL BE UP AND WALKING IMMEDIATELY AFTERWARD!

At the Vein and Laser Center of Northern Colorado, we pride ourselves on doing a good job for our patients. Schedule a free screening with us today, and we can see if EVLT is the right varicose veins treatment for you.



# Radial fiber ile hemoroid lazer tedavisi



# Applications

Colemed  
pharma



Colposcopy Applications  
HPV lesions  
Cysts of the mucosa  
Cervical intraepithelial neoplasia

## Vaginal Tightening Function

### Technical Data:

- Designed for performing minimally invasive gynecological treatments
- Less invasive, More effective
- Wavelength: 808nm/980nm/1470nm
- Power: 30W
- Length: 219.5mm
- Outside diameter: 35mm
- Inside diameter: 24mm






PLDD Solution for  
Neck & Beck Pain



Colemed  
pharma

Perkütan Lazer Disk Dekompresyonu  
PLDD lazer tedavisi





**EAR**



**NOSE**



**THROAT**


# Kulak Burun Boğaz KBB lazer tedavisi



## Application



Ear  
Cysts  
Accessory auricle  
Tympanitis  
Hemangioma



Nose  
Papilloma  
Epistaxis  
Stenosis & synechia  
Sinus surgery



Throat  
Uvulopalatoplasty  
Vocal cord polyps  
Strictures  
Epiglottectomy

# 4 tedavi başlıklı diş tutacağı lazer tedavisi

Başvuru:

Yumuşak doku sorunları, İmplantın iltihap önleyici etkisi,  
Kanal tedavisi, Periodontal cep tedavisi



Bio-stimulation



Therapy (TMJ)



## Application



Soft tissue cutting



Teeth whiteing



## 6. Profesyonel çoklu ekran işletim sistemi

- a) Sistemde seçenek olarak 10 adet fonksiyonumuz bulunmaktadır. İstedığınız zaman yeni bir el aleti olarak başka fonksiyonlar ekleyebilirsiniz.
- b) Tedavi süresini ayarlayabileceğiniz Zamanlayıcı vardır.
- c) Fizik tedavi fonksiyonunda hafıza fonksiyonu, sisteme 30 müşteri verisi kaydedebilirsiniz.

**Colemed  
pharma**

### Diode Laser System

Cutting Mode

EVLT Mode

ENT Mode

Dental Mode

Intimate Mode

PLDD Mode

### 980nm+1470nm Diode laser

1	2	3	4	5	6	7	8	9	10
11	12	13	14	15	16	17	18	19	20
21	22	23	24	25	26	27	28	29	30

### Diode Laser System

Vascular Removal

Anti-cellulite

Physiotherapy

Onychomycosis



OTHER  
FUNCTIONS

### 810nm Laser

### ARFURLA LASER SYSTEM

#### Physiotherapy

Energy

W

Pulse

ms

Frequency

Hz

Timer

CW Pulse Mode

Pulse Mode

Single

J-

min

Standby

## **Tırnak mantarı giderme, fizik tedavi ve damar temizleme lazerinizin avantajları:**

1. 810nm, 980nm, 1064nm,1470nm veya 980nm+1470nm gibi seçenek için çoklu dalga boyları
2. Yüksek güç, Sınıf IV
3. Almanya lazer çıkışının istikrarlı ve uzun ömürlü olmasını sağlamak için Lazer diyotu ithal etti
4. Almanya'dan ithal edilen el aletinin yan tarafındaki lens, optik iletim sırasında enerji kaybını azaltır
5. Kablo ve başlık konektörü arasında boşluk yok, fiberi tozdan koruyun, fiberi daha uzun ömürlü hale getirin
6. İkiden fazla fonksiyona ihtiyacınız varsa, her el aletinde Ayrı kablo bulunur, ekstra kablo için ödeme yapmanız gerekmez. Böylece farklı işlemlere ihtiyaç duyduğunuzda tedavi başlığını değiştirmenize gerek kalmaz.
7. Fizik tedavi için el aleti, El aleti üzerindeki ayak anahtarı ve kontrol düğmesi, bir lazer üzerinde iki kontrol sistemi ve 20mm-30mm arasında ayarlanabilir nokta boyutu.
8. Vasküler çıkarma el aleti için, seçenek için 8 nokta boyutu vardır, tedavide ihtiyacınıza göre seçebilirsiniz. Diğerleri piyasada sadece 3-4 spot büyüklüğüne sahiptir.
9. İşlev sayfalarında seçenek için 10 işlev
10. Optik fiber, seçenek için 3m 400um/600um'dur. Radyal fiber, keskin uçlu ve yuvarlak uçlu 3m, 600um'dur

# Teşekkürler

Hoş işbirliğimizi bekleyin!